

CATALOGUE 2024-25





5 Point - Calibration for pH & Conductivity)

Mode
pH
ORP
Conductivity
TDS
Salinity
Resistivity

°C





Water Quality Multiparameter LMMP30

pH MODE: Up to 5 points push-button calibration with auto-buffer recognition. Selectable pH buffer standards, (USA, DIN, NIST, User-defined) resolution (0.1, 0.01, 0.001pH)temperature units (°C or °F) Automatically calculates and displays electrode slope that helps you decide whether to replace pH electrode. Calibration due reminder prompts you to calibrate the meter regularly.

ORP MODE: Absolute or Relative mV measurements provides accurate ORP measured value .1 point offset calibration allows user adjusts displayed value to a known standard.

CONDUCTIVITY/TDS/SALINITY/RESISTIVITY MODE: Up to 5-point push-button calibration allows user to uses customized calibration solutions. Default 10uS, 84uS, 1413uS, 12.88mS and 111.8mS)Adjustable TDS conversion factors, cell constants, temperature coefficients and normalization temperature. (20°C or 25°C)Allows user selects linear or pure water temperature compensation modes. In salinity measurement modes, selectable practical salinity (PSU) or seawater measurement standards. (ppt)

SALIENT FEATURES

- Selectable either temperature units between degrees Celsius (°C) and degrees Fahrenheit (°F)
- Automatic Temperature Compensation (ATC) ensures high accuracy of measurements
- Automatic endpoint function freezes the stable measured values for easy reading and recording data.
- ♣ Help message as a operational guide to helps you understand how to use meter.
- Setup menu lets user customize the calibration points. stability conditions of measurement. temp. units and date/time to meet personal preference.
- Reset function allows user to returns all setting to the factory default values.
- Expanded memory stores and recalls upto 500 readings.
- Built-in real-time clock stamps stored data to meets GLP standards.

STANDARD ACCESSORIES: Main Instrument, pH Electrode, pH Buffer - (4.01, 7.00pH), Conductivity Cell 1K, Calibration Solution, Electrode Stand, ATC, USB Adaptor, Test Certificate







CATALOGUE 2024-25

Multiparameter	HSN Code : 90275090
Model	LMMP30
pH Range	-2.000~20.000pH
Resolution	0.001, 0.01, 0.1pH, Selectable
Accuracy	±0.002pH
Calibration Points	Up to 5 points
pH Buffer Option	USA, (pH1.68, 4.01,7.00,10.01,12.45), NIST (pH1.68/4.01/6.86/9.18/12.45), DIN (pH1.09,4.65,6.79,9.23,12.75)
ORP (mV) Range / Calibr	0 ~ ±1999.9mV / 1 point
Resolution	Q.1, /1mV, Selectable
Accuracy	±0.2mV
Measuring Mode	Absolute mV or Relative mV
Conductivity Range	0~200mS/cm
Resolution	0.01 / 0.1/ 1
Accuracy	±0.5% Full Scale
Cell Constant	0.1/1/10/User defined 0.0~10.0% per °C
Temperature Coefficient	0.0~10.0% per °C
Calibration Points	Up to 5 points
Calibration Solutions	10uS/cm, 84uS/cm, 1413uS/cm, 12.88mS/cm, 111.8mS/cm
Compensation Modes	Linear or Pure Water
TDS Range	0~10ppt (Max.200 ppt, Depending on factor setting)
Resolution	0.01/0.1/1
Compensation Modes TDS Range Resolution Accuracy TDS Factor	±1% Full Scale
TDS Factor	0.0 ~1.0
Salinity Range	0~80ppt
Resolution	0.01 / 0.1 /1
Accuracy	±1% Full Scale
Measurement Mode	Practical Salinity (PSU) or Seawater Measurement Standard (ppt)
Resistivity Range	0~100ΜΩ
Resolution	1/0.1/0.01
Accuracy	±1% Full Scale
Temp. Range	0~100°C, / 32~212°F
Temp. Compensation	0~100°C, / 32~212°F Automatic or Manual
Temp. Accuracy	±0.5°C, ±0.9°F
Auto-Off Function	10 / 20 / 30 minutes after last key pressed
Reset & Auto Hold Functi	n Yes
Electrode Diagnosis	Slope, Offset
Stability Conditions	Slope, Offset High or Low, Selectable 1 to 30 days, Selectable
Calibration due reminder	1 to 30 days, Selectable
Memory	500 data sets with Date/Time
Connector	(pH) BNC 6-pin connector - (Conductivity)
Dimensions (LxWxH)& W	ight 210 × 188 × 60 mm & 1.2kg