





**Automatic Elisa Plate Washer
AEPW99**



SALIENT FEATURES

- ✦ 4.3" colour Touch screen
- ✦ Process Mode - Full plate wash / row by row wash
- ✦ Plate washing effect is excellent and highly efficient.
- ✦ Residual liquid cleaning function
- ✦ Compatible with various microplate types
- ✦ Flat & bottom U-shaped, V-shaped, and C-shaped wells
- ✦ Compact Design, High cleaning efficiency
- ✦ Liquid pump type - Diaphragm pump

Automatic Elisa Plate Washer AEPW99

Automatic ELISA Plate Washer streamlines the ELISA process by efficiently washing microplates, reducing Manpower and minimizing errors. It ensures precise liquid handling, uniform washing, and contamination control, making it an essential tool for high-throughput laboratories.



E - Catalogue

HSN Code : 90275090

SPECIFICATIONS

MODEL	AEPW99
Make	LABMAN
Residual liquid volume	$\leq 1.5\mu\text{L}$
Plate Hole	96
Wash head	8 Channels/ 12 Channels cleaning head
Number of washes	1-99
Washing Volume	50ml - 3000ml
Wash bottles	3 bottles
Injection accuracy	$\leq 2\%$ (300ul)
Soaking time	1S - 99min
Timer	Available
Display	Touch Screen
Packing Dimensions(L x W x H)	453 x 388 x 155 mm
Gross Weight	14Kg

★ Specifications are subject to change without notice.

STANDARD ACCESSORIES : Main Instrument, Washing Head, Wash Bottles, , Power cord, Instruction Manual, Test Certificate

LABMAN Product Range 2025-26

LMFD61, 62, 63,64  Lyophilizer Freeze Dryer (-60°C)	LMFD83, 84  Lyophilizer Freeze Dryer (-80°C)	LMRE25  Rotary Evaporator	ADM3  Automatic Density Meter	MSOD6  Microvolume Spectrophotometer
NAE32  Nucleic Acid Extractor	EMR4  Elisa Microplate Reader	GPTC96  Gradient PCR Thermal Cycler	HO30  Hybridization Oven	UVCL254  UV - Crosslinker
TGH12, 24  Tissue Grinder Homogenizer	PBH80T  Paddle Blender Homogenizer	TOC2NDR  Total Organic Carbon Analyzer	TOC1CD  Total Organic Carbon Analyzer	LMDS600, 1200  Microwave Digestion System
HPLC1IS, HPLC2GS  High Performance Liquid Chromatography	APT1800  Automatic Potential Titrator	KFT250VC  Karl Fischer Titrator	COD1500  COD Analyzer with Digestor	BOD4000  BOD Analyzer
MTVM50  Multi Tube Vortex Mixer	LMVM60D  Vortex Mixer	OS20  Orbital Shaker	TRM6  Tube Roller Mixer	DBNE12, 24  Dry Bath Nitrogen Evaporator





LABMAN Product Range 2025-26

DBI11T  Dry Bath Incubator	DBI15SH, DBI16SC  Dry Bath Incubator	IRS14, 35  Infrared Sterilizer	GGA200  Air Generator	GGA301, 302  Hydrogen Generator
GGN303  Nitrogen Generator	GHA304  Hydrogen Air Generator(2 in 1)	LMVP10  Vacuum Pump	LMVP20  Vacuum Pump	LMCO1  Column Oven
LMIO6  Mini Thermostatic Incubator	GHP400  Graphite Hotplate	LOBMS5  Oil Bath Magnetic Stirrer	LMMS4LC, 20LC  Ceramic Hotplate Magnetic Stirrer	LMMS300  Hotplate Magnetic Stirrer
MMS15S  Multiposition Magnetic Stirrer	MMS12H  Multiposition Magnetic Stirrer	WBMS3  Water Bath Magnetic Stirrer	LMWB2H, 4H, 6H, 8H  Digital Water Bath	HOS20L  Overhead Stirrer
LMHH160  Hand Held Homogenizer	PS156S  Peltier / Sipper System	LUV6HG  PMT Spectrophotometer	LUV2000T, 2100T  Single / Double Beam UV-Vis. Spectrophotometer with 21CFR	LMSPUV SERIES  Double Beam UV-Vis. Spectrophotometer

LABMAN Product Range 2025-26

LMSPUV1200  UV-Vis. Single Beam Spectrophotometer	LMSPV320, 325  Visible Spectrophotometer	ADP100CFR, 95CFR  Automatic Digital Polarimeter with CFR	ADP45, 50  Automatic Digital Polarimeter (±45°)	ADP75, 90  Automatic Digital Polarimeter (±90°)
ADP40  Automatic Digital Polarimeter	LMPM180  Manual Polarimeter	MAR100  Manual ABBE Refractometer	RFM100CFR  Digital Refractometer CFR	RFM950, 970  Digital Refractometer
LMAR1317T  ABBE Refractometer	LMDV300H  Digital Viscometer with Thermoset	LMDV60, 100, 200  Digital Viscometer (L Series)	LMDV400R  Digital Viscometer (R Series)	LMKV2  Krebs Viscometer
LMDV111, 222  Digital Viscometer with SSA	LVA30  Low Viscosity Adaptor	LMCR80F  Visual Colorimeter	LMUC2D  Ultrasonic Cleaner	LMUC3D, 6D, 9D, 12D, 25D  Ultrasonic Cleaner - Degas
LMUC50, 100  Industrial Ultrasonic Cleaner	LMUC13M, 35M  SS304 - Ultrasonic Cleaner	LMPH15  pH Meter with Printer	LMPH12, LMPH10H  pH Meter	LMPH9  pH Meter - Basic

LABMAN Product Range 2025-26

LMCM20H  Conductivity / TDS Meter	LMMP30, 32  Multiparameter (7 in 1 / 9 in 1)	LMION40  Ion Meter	LMDO50  Dissolved Oxygen Meter	LMHD60  Water Hardness Meter
OPH11  Online pH Controller	OCM21  Online Conductivity Controller	LMTB200  Turbidity Meter	RCB1010T, 1020T, 1030T  Refrigerated Circulating Bath	RCB620T  Refrigerated Circulating Bath
RWC320  Recirculating Water Chiller	LMDO4P  Portable DO Meter	LMPH1P  Portable pH / ORP Meter	LMCM2P  Portable Conductivity / TDS Meter	LMTB5P  Portable Turbidity Meter
PRO656, 900, 1200, 1800  Probe Sonicator	PRO660R  Probe Sonicator with Remote	PRO650  Probe Sonicator	LMIF30T, 50T, 100T  Ice Flaker	LMIF200T  Ice Flaker
LMIM15  Ice Maker	LMIM20C  Ice Maker with Crusher	ISE19, 62  Nitrate / Fluoride ISEs	ER SERIES  pH Electrode / DO Electrode	CC1, CC2, CC3  Conductivity Cell