



Catalogue

GGA200

Air Generator



13 MONTHS
Guarantee

2024-25

GGA200

CATALOGUE 2024-25

**Air Generator****GGA200****SALIENT FEATURES**

- ✦ Air generator is to provide pure compressed air source device
- ✦ Generator using phase process control technology and gas condensation separation purification technology step by step.
- ✦ Dual pressure stabilization
- ✦ Stainless Steel Gas Storage Tank
- ✦ Silicon Rubber Ring with Low Sulphur Content
- ✦ High Pressure Resistance, High Precision Stabilization
- ✦ Adsorbent Molecular Sieves regeneration
- ✦ Color changing silicon should be heated at 120-140°C until it turns completely blue.
- ✦ High purity and compressed air flow is very stable

HSN Code : 84051010



E - Catalogue

SPECIFICATIONS	MODEL	GGA200
	Make	LABMAN
	Maximum Flow	2000mL/min
	Output Gas Pressure	0.4MPa
	Color changing Silicon	120-140°C
	Output Flow Display	LED Display
	Phase	1 Phase
	Power Wattage	180W
	Power supply	AC 220V / 50Hz
	Dimension (LxWxH) & Weight	565 x 400 x 500 mm
	Weight	29.8kg

★ Specifications are subject to change without notice.

STANDARD ACCESSORIES : Main Instrument, O-Ring, Rubber Bush, Air Intake Screw, Power Cord, Instruction Manual, Test Certificate

LABMAN PRODUCTS

Laboratory & Scientific Instruments



				
Total Organic Carbon Analyzer	Microwave Digestion System	COD Analyzer with Digestor	BOD Analyzer with CFR	Bio Chemistry Analyzer
				
Touch Screen Digital Polarimeter with CFR	Touch Screen Digital Refractometer with CFR	Touch Screen ABBE Refractometer	Double / Single Beam UV-Vis.Spectrophotometer	Double Beam UV-Vis.Spectrophotometer
				
Single Beam UV-Vis.Spectrophotometer	Visible Spectrophotometer	Peltier / Sipper System	Gradient Thermal Cycler PCR	Visual Colorimeter
				
Ice Flaker - Touchscreen	Refrigerated Circulating Bath	Refrigerated Circulating Bath	Industrial Ultrasonic Cleaner	Degas - Function Ultrasonic Cleaner
				
Probe Sonicator with Remote	Ceramic Hotplate Magnetic Stirrer	Hotplate Magnetic Stirrer	Digital Viscometer	Low Viscosity Adaptor
				
Dry Bath Incubator	Vortex Mixer	Digital Water Bath	pH meter with Printer	pH/Cond./Multipara./Ion/DO Meter
				
Turbidity Meter	Online pH/Conductivity Controller	Nitrate / Fluoride - ISE	Conductivity Cell	Probe for Probe Sonicator