



Catalogue

GHA304

Hydrogen Air Generator (2 in 1)



13 MONTHS
Guarantee

2024-25

GHA304

CATALOGUE 2024-25



SALIENT FEATURES

- Hydrogen Air Generator is a multi-purpose instrument
- LED Display, Easy to Operate
- High degree of automation, gas can be produced by starting the switch.
- Over pressure protection with anti alkali liquid is used
- Two-stage overvoltage protection & voltage stabilization devices
- Drainage switches & dust filtration for multi-stage purification
- Low Vibration, Low Noise, Stable output pressure & flow rate.

Hydrogen Air Generator (2 in 1)

GHA304

HSN Code : 84051010



E - Catalogue

SPECIFICATIONS	MODEL	GHA304
	Make	LABMAN
	Air quality	Dry, oil-free
	Maximum H ₂ Flow	300ml/min.
	Maximum Air Flow	2000ml/min.
	Gas purity	99.999%
	Output H ₂ Gas Pressure	0.3mpa
	Output Air Gas Pressure	0.4mpa
	Display	LED Display
	Environmental temperature	10°C -40°C
	Power Wattage	330W
	Power supply	AC 220V / 50Hz
	Dimension (LxWxH) & Weight	525 x 435 x 750 mm
	Weight	41kg

★ Specifications are subject to change without notice.

STANDARD ACCESSORIES : Main Instrument, O-Ring, Rubber brush, Potassium Hydroxide(KOH)120g, Air Intake Screw, Power Cord, Instruction Manual, Test Certificate.

LABMAN PRODUCTS

Laboratory & Scientific Instruments



				
Total Organic Carbon Analyzer	Microwave Digestion System	COD Analyzer with Digestor	BOD Analyzer with CFR	Bio Chemistry Analyzer
				
Touch Screen Digital Polarimeter with CFR	Touch Screen Digital Refractometer with CFR	Touch Screen ABBE Refractometer	Double / Single Beam UV-Vis.Spectrophotometer	Double Beam UV-Vis.Spectrophotometer
				
Single Beam UV-Vis.Spectrophotometer	Visible Spectrophotometer	Peltier / Sipper System	Gradient Thermal Cycler PCR	Visual Colorimeter
				
Ice Flaker - Touchscreen	Refrigerated Circulating Bath	Refrigerated Circulating Bath	Industrial Ultrasonic Cleaner	Degas - Function Ultrasonic Cleaner
				
Probe Sonicator with Remote	Ceramic Hotplate Magnetic Stirrer	Hotplate Magnetic Stirrer	Digital Viscometer	Low Viscosity Adaptor
				
Dry Bath Incubator	Vortex Mixer	Digital Water Bath	pH meter with Printer	pH/Cond./Multipara./Ion/DO Meter
				
Turbidity Meter	Online pH/Conductivity Controller	Nitrate / Fluoride - ISE	Conductivity Cell	Probe for Probe Sonicator