

LMDV60, 100, 200

CATALOGUE 2024-25



Digital Rotational Viscometer

LMDV60,100,200

HSN Code : 90278910

SALIENT FEATURES

- ✦ Auto Range function
- ✦ Selectable Speed
- ✦ Direct Viscosity reading and temperature display
- ✦ Big Ultra Bright Backlight LCD Display

Digital Direct Reading Rotational Viscometer

Labman LMDV series is Rotational Digital Direct Reading Viscometer to measure absolute viscosity of Newton Liquids as well apparent viscosity of non – Newton liquid featured by high flexibility reliable Test result, easy operation and good appearance

Optional Accessories :

- Low Viscosity Adaptor
- Small Sample Adaptor
- Refrigerated Circulating Bath



E - Catalogue

SPECIFICATIONS

Model	LMDV60	LMDV100	LMDV200
Make	LABMAN		
Measuring Range mpa.s	20-2,000,000	20-2,000,000	20-6,000,000
Speed	0.3 to 60	0.3 to 100	0.1 to 200
No of Speed	8	18	54
No of Standard Spindle	4		
Display	LCD		
Temperature Monitoring	Yes		
Measurement Accuracy	+ / - 1% (FS)		
Repeatability	+ / - 0.5 % (FS)		
Output	RS232C Interface		
Power Supply	Adaptor, Using AC 220V,50/60Hz		
Dimension LxWxH (mm)	400 x 350 x 360		
Weight	9 kg		10 kg

★ Specifications are subject to change without notice.

STANDARD ACCESSORIES : Main Instrument, Protection Cap, DC Adaptor, Spindle Set - L1, L2, L3, L4, RTD Probe, Test Certificate, Spindle Guard & Wrench, Indicator Stand (LMDV60)