

LMIF200T

CATALOGUE 2024-25



Ice Flaker - Touchscreen
LMIF200T

SALIENT FEATURES

- ◆ Touchscreen Display, Timer Function.
- ◆ Storage Bin Detachable
- ◆ Fully Automatic microprocessor control
- ◆ Noiseless operation, Temperature Function.
- ◆ Compact Design and low maintenance.
- ◆ Continuous Ice Flakes Output.
- ◆ CFC free compressor, Base feet 8 nos.
- ◆ Stainless Steel construction.
- ◆ Low water level Indication.
- ◆ Front & back Ventilation Design.

APPLICATION : Hospitals, Educational Research
Institutes & Biotechnology Laboratories



E - Catalogue

HSN Code : 84186990

MODEL	LMIF200T
Make	LABMAN
Capacity (kg per 24hours)	200 kg
Storage Capacity	60kg
Display Type	Touch screen
Ice Type	Flake ice
Storage Bin Material	ABS
Cooling Method	Compressor cooling with fan cooling
Ambient Temperature	Below 25°C for Optimum Working
Water Temperature (Recommended)	20°C
Construction Material	Stainless steel
Noise Level & Water Supply Size	58db & 25mm
Refrigerant	R290A
Auto off	yes
Input Power	850W
Packing Dimension (LxWxH)	790 x 715 x1479 mm
Net Weight	68kg
Gross weight	78kg

STANDARD ACCESSORIES : Main Instrument, Water Filter with Inlet & Outlet Tube, Inlet Tube, Drainage Tube, Ice Shovel (Scoop), Leg, Test Certificate.

★ Specifications are subject to change without notice.