



Catalogue

LMIO6

Mini Thermostatic Incubator



2024-25

LMIO6

CATALOGUE 2024-25



Mini Thermostatic Incubator

LMIO6

Applications: It is used for industrial and small enterprises, health, Biochemistry, agricultural science and other Scientific Research Departments to do constant temperature Culture experiments such as bacteria Culture, Breeding, Fermentation etc.

HSN Code : 90321090



E - Catalogue

SPECIFICATIONS	MODEL	LMIO6
	Make	LABMAN
	Volume	6L
	Power	100W
	Temp.range	RT +5 ~ 45°C
	Ambient temperature	5-30°C
	Temp. Resolution & Accuracy	0.1°C & 1°C
	Shape size	100 x 125mm
	Structure	Portable
	Time Setting	Yes
	Inner chamber size (L×W×H)	215 x 220 x 205 mm
	Outer size mm (L×W×H)	300 x 320 x 355 mm
	Power supply	AC 220V / 50Hz
	Dimension (L×W×H)	400 x 390 x 420 mm
	Gross weight	10kg


★ Specifications are subject to change without notice.

STANDARD ACCESSORIES : Main Instrument, SS Tray, Power Adaptor, Instruction Manual, Test Certificate.

LABMAN PRODUCTS

Laboratory & Scientific Instruments



APT1800  Automatic Potential Titrator	KFT250VC  Karl Fisher Titrator	KFT1B  Karl Fisher Titrator - Basic	SPB80T  Stomacher Paddle Blender	GGA200  Air Generator
GGH301. 302  Hydrogen Generator	GGN303  Nitrogen Generator	GHA304  Hydrogen / Air Generator (2 in 1)	LMVP10  Vacuum Pump (Basic)	LMVP20  Vacuum Pump (Chemical Resistance)
LMCO1  Column Oven	LMIO6  Portable Thermostatic Incubator	GHP400  Graphite Hotplate	LOBMS5  Oil Bath Magnetic Stirrer	LMMS4LC, 20LC  Ceramic Hotplate Magnetic Stirrer
TOC1CD  Total Organic Carbon Analyzer	TOC2NDIR  Total Organic Carbon Analyzer	LMDS600, 1200  Microwave Digestion System	ADP100CFR  Automatic Digital Polarimeter CFR	ADP45, 50  Automatic Digital Polarimeter
ADP75, 90  Automatic Digital Polarimeter	ADP95CFR  High Speed Automatic polarimeter	RFM100CFR  Automatic Digital Refractometer	RFM950, 970  Digital Refractometer	LMAR1317T  ABBE Refractometer
PCR9602G  Gradient Thermal Cycler PCR	RCB1020T, 1030T  Refrigerated Circulating Bath	RCB620T  Refrigerated Circulating Bath	LMIF30T, 50T, 100T  Ice Flaker - Touch Screen	LMIF200T  Ice Flaker - Touch Screen
LUV2000T, LUV2100T  Double Beam UV-VIS Spectrophotometer	LMSPUV1900, 1908, 1920, 1928  Double Beam UV - VIS Spectrophotometer	LUV6HG  PMT Spectrophotometer	LMSPUV1200, 1000B  Single Beam UV - VIS Spectrophotometer	LMSPV325, 320  Visible Spectrophotometer

LABMAN PRODUCTS

Laboratory & Scientific Instruments



LMUC3D, 6D, 9D, 12D, 25D



Digital Ultrasonic Cleaner
PRO656, 900, 1200, 1800

LMUC35M



Digital Ultrasonic Cleaner
PRO660R

LMUC50, 100



Industrial Ultrasonic Cleaner

COD1500



COD Analyzer with Digestor

BOD4000



BOD Analyzer



Touch Screen
Probe Sonicator

DBI11T



Touch Screen
Probe Sonicator with Remote

DBI18, 20



Probe

LMDV300H



Ceramic Hotplate
Magnetic Stirrer

LMMS5LC



Hotplate
Magnetic Stirrer

LMMS300



Touch Screen
Dry Bath Incubator

DBI11T



Dry Bath Incubator

DBI18, 20




Viscometer with Thermoseal

LMDV300H



Digital Viscometer

LMDV60, 100, 200



Krebs Viscometer

LMKV3



pH Meter- Economic

LMMP32



pH Meter- 3 Points, 5 Points

LMION40



pH Meter With Printer

LMPH15



Conductivity Meter - 3 Points

LMCM20H



Multiparameter - 7 in 1

LMMP30



Multiparameter (9 in 1)

OCM21



Ion/ pH Meter

LMION40



Dissolved Oxygen Meter

LMD050



Water Hardness Meter

LMHD60



online pH Controller

OPH11



Online Conductivity Controller

OCM21



Turbidity Meter

LMTB200



Visual Colorimeter

LMCR80F



Digital Water Bath

LMWB04



Small Sample Adaptor
Low Viscosity Adaptor

LVA30, SSA15



pH Electrode, pH/ ORP Electrode

ER1-10



Conductivity Electrode

CC1, 2, 3



Fluoride

Nitrate

ISE19

ISE62



Vortex Mixer

LMVM60D



Peltier / Sipper System

PS1565