

**LMMS5LC**

CATALOGUE 2024-25



### SALIENT FEATURES

- ◆ Digital LCD display with back light
- ◆ Big Size Ceramic hot plate with flame protection and durability
- ◆ Heating and stirring can proceed simultaneously.
- ◆ Temperature & Stirring speed is stepless adjusted.
- ◆ Real temperature & speed display on the screen simultaneously
- ◆ Two sensors, outside PT 100 sensor controls the liquid temperature, inside K type sensor control the hot plate temperature (selectable)
- ◆ Over temperature alarm

## Ceramic Hotplate Magnetic Stirrer

**LMMS5LC**



E - Catalogue

HSN Code : 84798200

MODEL	LMMS5LC
Make	LABMAN
Display	Digital
Speed (RPM)	200-2000
Plate Size (mm)	190 x 190
Top Plate Type	Ceramic
Max. Volume (ml)	5000
Temp. Controller	PID
Max Temp. ( Plate Surface)	300°C
Temp.resolution	0.1°C
Temp. Accuracy	±0.5°C(<100°C)±1°C(<300°C)
Power Wattage	600 W
Power supply	AC 220V,50/60Hz
Dimension LxWxH (mm)	300 x 190 x 126
Weight	5 kg

**SPECIFICATIONS**

★ Specifications are subject to change without notice.

**STANDARD ACCESSORIES** : Main Instrument, Teflon Paddle(Stirrer), ATC, Stand Set, Power Cable, Test Certificate