

**LMPH9**

CATALOGUE 2024-25



**pH Meter - Economic**  
**LMPH9**

## SALIENT FEATURES

- Manual Temperature Compensation
- There is a microprocessor chip inside the meter and the appearance of meter looks nice with user friendly operation. Some of the features are as under
- Large LCD display with blue background simultaneous display of pH and temperature
- It has intelligent functions for calibration, data storage / Memory (upto 25 prints)
- Automatically display the percentage of slope after calibration is over.
- Equipped with electrode holder and calibration buffer solution, which is more convenient to use.



E - Catalogue

HSN Code : 90275090

MODEL	LMPH9
Make	LABMAN
Range	0 to 14.00 pH
pH Resolution	0.1 pH / 0.01 pH
Accuracy	± 0.01 pH
Calibration Points	up to 3 points
pH Buffer options	USA (pH 4.01, 7.00, 10.01)
mV Range	±1999 mV
Resolution & Accuracy	± 1 % FS ± 1 digit & 1mV
Temperature Range	0°C to 100°C
Accuracy	± 0.5°C
Resolution	0.1°C
Temperature Compensation	Manual
Connector	BNC
Display	LCD
Memory	Upto 25 prints
Power Requirements	Adaptor, Using AC 220V,50/60Hz
Packing Dimensions ( Lx W x H)	160 × 190 × 170 mm
Weight	750g

SPECIFICATIONS

★ Specifications are subject to change without notice.

**STANDARD ACCESSORIES :** pH Meter, pH Electrode, pH Buffer - 4.01, 7.00, Power Adaptor, Electrode Stand, Test Certificate.