

**LMVM60D**

CATALOGUE 2024-25



**Vortex Mixer  
LMVM60D**

**SALIENT FEATURES**

- LCD Display with white back light
- Speed Timer (1~99h 59m 59s)
- Operating Mode - 2 Types - Auto run sensor, continuous run used for various Mixing application
- Supply with Standard Head
- Soft silicone foot pads
- Shockproof, Brushless Motor
- Aluminium Alloy Base stable body

Digital display vortex mixer is powered by a brushless motor and features an eccentric structural design.

This innovative mixer effectively blends and oscillates liquids within a container. It finds extensive applications in various fields, including biochemistry, genetic engineering, and medicine for conducting experiments

HSN Code : 84798200



E - Catalogue

Model	LMVM60D
Make	LABMAN
Speed	100-2800rpm (adjustable)
Display	LCD with white back light
Display Function	Speed & Time
Adjustment method	Keypad
Operation Mode	2 Types - Auto run sensor, continuous run
Oscillation mode	Circumference
Base material	Aluminium alloy
Shaking Movement	Orbital
Orbital diameter	4mm
Ambient Temperature	5°C-35°C
Humidity	≤ 80%
Input power	60W
Power supply	AC 220V,50/60Hz
Dimension (LxWxH) & Gross weight	250x220x190mm & 3.5kg

**SPECIFICATIONS**

**STANDARD ACCESSORIES** : Main Instrument with Standard Head, Power Cable with Adaptor, Test Certificate

★ Specifications are subject to change without notice.