



Catalogue

LOBMS5

Oil Bath Magnetic Stirrer



2 YEARS
Guarantee

2024-25

LOBMS5

CATALOGUE 2024-25



SALIENT FEATURES

- ◆ 5L Volume Capacity, PID temperature control Technology
- ◆ LED Display, Stepless Speed regulation
- ◆ Low Noise, Stable Operation
- ◆ PT100 Temperature Sensor
- ◆ High Temperature Resistance
- ◆ Fast Heating Speed
- ◆ Uniform temperature and High efficiency
- ◆ Strong Magnet to ensure Sufficient Suction
- ◆ Heating Method - Electric & Tube Heating

Oil Bath Magnetic Stirrer

LOBMS5

HSN Code : 84798200



E - Catalogue

MODEL	LOBMS5
Make	LABMAN
Maximum Volume	5000ml
Display	LED
Speed range	0-2800 RPM (Adjustable)
Stirring Power	12 W
Liner Height	145mm
Liner Material	SS304
Temperature Sensor, Control	PT100, PID Temperature
Time Setting	Yes
Max. Temp,	300°C
Temp. Resolution	±1°C
Power Wattage	1000 W
Power supply	AC 220V / 50Hz
Dimension (LxWxH) & Weight	290 x 290 x 260 mm & 6kg

SPECIFICATIONS

★ Specifications are subject to change without notice.

STANDARD ACCESSORIES : Main Instrument, Pot Lid 5L, Stir Bar, Temperature Sensor, Column Stand, Temperature Sensor Bracket, Clamp, Beaker Holder, Cross Clamp, Instruction Manual, Test Certificate.

LABMAN PRODUCTS

Laboratory & Scientific Instruments



				
Total Organic Carbon Analyzer	Microwave Digestion System	COD Analyzer with Digester	BOD Analyzer with CFR	Bio Chemistry Analyzer
				
Touch Screen Digital Polarimeter with CFR	Touch Screen Digital Refractometer with CFR	Touch Screen ABBE Refractometer	Double / Single Beam UV-Vis.Spectrophotometer	Double Beam UV-Vis.Spectrophotometer
				
Single Beam UV-Vis.Spectrophotometer	Visible Spectrophotometer	Peltier / Sipper System	Gradient Thermal Cycler PCR	Visual Colorimeter
				
Ice Flaker - Touchscreen	Refrigerated Circulating Bath	Refrigerated Circulating Bath	Industrial Ultrasonic Cleaner	Degas - Function Ultrasonic Cleaner
				
Probe Sonicator with Remote	Ceramic Hotplate Magnetic Stirrer	Hotplate Magnetic Stirrer	Digital Viscometer	Low Viscosity Adaptor
				
Dry Bath Incubator	Vortex Mixer	Digital Water Bath	pH meter with Printer	pH/Cond./Multipara./Ion/DO Meter 2mm to 25mm
				
Turbidity Meter	Online pH/Conductivity Controller	Nitrate / Fluoride - ISE	Conductivity Cell	Probe for Probe Sonicator